

SARS-CoV-2 and Influenza A+B Antigen Rapid Test

FOR PROFESSIONAL USE ONLY

Product Name

SARS-CoV-2 and Influenza A+B Antigen Rapid Test

Intended Use

This product is used for in vitro qualitative detection of SARS-CoV-2 and influenza A virus, influenza B virus antigen in human nasopharyngeal swab or oropharyngeal swab samples.

This reagent is only used in clinical laboratory, medical institutions or real-time inspection by professional medical personnel, not suitable for family test, the test results are only for clinical reference, it is recommended to conduct comprehensive analysis of the disease condition in combination with clinical manifestations of patients and other laboratory tests; it is not suitable for screening of general population.

Test Principle

According to the gold immunochromatographic test principle, double antibody sandwich method was used to detect SARS-CoV-2 antigen and influenza A virus antigen, influenza B virus antigen in the samples.

When the sample contains SARS-CoV-2 antigen, it forms a complex with the gold-labeled antibody (SARS-CoV-2 monoclonal antibody 1). The complex moves forward under the action of the chromatogram and combines with the coating antibody (SARS-CoV-2 monoclonal antibody 2) at the T line to form a complex and develop color (T line), the result is positive. When there is no SARS-CoV-2 antigen in the sample, the T line does not form a complex, and there is no red band, which is a negative result.

When the sample contains influenza A virus antigen, it forms a complex with the gold-labeled antibody (influenza A virus monoclonal antibody 1). The complex moves forward under the action of the chromatogram and combines with the coating antibody (influenza A virus monoclonal antibody 2) at the T line to form a complex and develop color (T line), the result is positive. When there is no influenza A virus antigen in the sample, the T line does not form a complex, and there is no red band, which is a negative result.

When the sample contains influenza B virus antigen, it forms a complex with the gold-labeled antibody (influenza B virus monoclonal antibody 1). The complex moves forward under the action of the chromatogram and combines with the coating antibody (influenza B virus monoclonal antibody 2) at the T line to form a complex and develop color (T line), the result is positive. When there is no influenza B virus antigen in the sample, the T line does not form a complex, and there is no red band, which is a negative result.

Regardless of whether the sample contains SARS-CoV-2 /influenza A virus/influenza B virus antigen, the gold-labeled biotinylated BSA will bind to the coated antibody at line C to form a complex and develop color (line C).

Main component

Detection card: the detection line is coated with antibodies against SARS-CoV-2 monoclonal antibody 2, influenza A virus monoclonal antibody 2, influenza B virus monoclonal antibody 2, Gold standard pad solid phase SARS-CoV-2 monoclonal antibody 1, influenza A virus monoclonal antibody 1, influenza B virus monoclonal antibody 1, Biotinylated BSA. The quality control line is coated with streptavidin-conjugated IgG.

Sample extract: Tris(hydroxymethyl)methyl aminomethane buffer with surfactant.

Swab and sample extraction tube are optional.

MATERIAL NEEDED BUT NOT PROVIDED

1. Timer
2. Personal protective equipment, such as protective gloves, medical mask, goggles and lab coat.
3. Appropriate biohazard waste container and disinfectants.

Storage And Shelf-Life

Store as packaged in the sealed pouch at 4-30°C, avoid hot and sunshine, dry place, valid provisional for 12 months. DO NOT FREEZE. Some protective measures should be taken in hot summer and cold winter to avoid high temperature or freeze-thaw. Do not open the inner packaging until ready, it must be used in one hour if opened (Humidity ≤ 60%, Temp: 20°C-30°C). Please use immediately when the humidity > 60%.

Sample Requirement

Sample Collection

Nasopharyngeal swab collection method:

The operator holds the swab by the right hand and holds the head of the subject fixedly by left hand. Do not overexert to avoid traumatic hemorrhage. When the cusp of the swab touching the paries posterior of the pharyngonasal cavity, letting the swab remain in the place for a few seconds (about 3 seconds) and rotating the swab gently for one cycle, and then remove the swab slowly. Using the same swab, repeat this process for the other nostril to ensure that an adequate sample is collected from both nasal cavities.

Collection method of oropharyngeal swab:

The head of the person to be collected is slightly tilted and his mouth is wide open, exposing the pharyngeal tonsils on both sides. Wipe the swab across the root of the tongue. Wipe the pharyngeal tonsils on both sides of the person to be collected back and forth with a little force for at least 3 times, and then wipe up and down the posterior pharyngeal wall for at least 3 times. Avoid touching your tongue, cheeks or teeth when sampling. Just after drinking water or beverages, sampling samples cannot be used for testing.


Note: The sample should not be inactivated.

Sample preservation

The samples of human nasopharyngeal swabs or oropharyngeal swabs should be placed in the sample extract immediately after collection and tested as soon as possible within 1 hour. Long term storage is not recommended.

Sample Treatment

Add 500µl (~20 drops) of sample extract to the 0.5 mark of the sampling tube, dip the swab after collecting the sample into the sample extract, make the sample extract fully permeate the swab, rotate and squeeze the swab 10 times, then pull out the swab, and take the stranded liquid as the sample to be tested.




Test Procedure

Instructions must be read entirely before taking the test. Leave the reagent and

sample at room temperature for 30 minutes before use. Return to room temperature. Do not open the inner packing until it is ready. Use it as soon as possible after opening the inner packing.

1. Open the tear hole of the aluminum foil bag, take out the test card and lay it flat.
2. Apply 2-3 drops of the treated sample extract (60µl-80µl) vertically into the sample well of the test cassette.
3. The results are observed after 15 minutes and showed no clinical significance after 20 minutes.




Interpretation Of Result

POSITIVE: Two distinct colored lines appear. One line should be in the control region (C) and the other line should be in the test region (T).

NEGATIVE: One colored line appears in the control region (C). No apparent colored line appears in the test region (T). The negative result does not indicate the absence of analytes in the sample, it only indicates the level of tested analytes in the sample is less than cut-off level.

INVALID: No colored lines in the position of control line appear, indicating that the operator error or reagent failure. Verify the test procedure and repeat the test with a new testing device.



Limitation

1. The result of the product should not be taken as a confirmed diagnosis, for clinical reference only. Judgement should be made along with RT-PCR results, clinical symptoms, epidemic condition and further clinical data.
2. If the virus antigen level in the sample is lower than the detection limit, the test result may be negative.
3. As the duration of the disease increases, the number of antigens in the sample may decrease. After the sample is collected, compared with RT-PCR analysis, five days after the onset of symptoms, the result may be negative.
4. Due to the limitation of the detection method, the negative result cannot exclude the possibility of infection. The positive result should not be taken as a confirmed diagnosis. Judgement should be made along with clinical symptoms and further diagnosis methods.
5. This reagent can only qualitatively detect SARS-CoV-2 antigen and influenza A virus antigen, influenza B virus antigen in human nasopharyngeal swab,

oropharyngeal swab. It cannot determine the certain antigen content in the samples.

6. The accuracy of the test depends on the sample collection process. Improper sample collection, improper sample transportation and storage or freezing and thawing of the sample will affect the test results.

7. It is optimum when eluting swabs with the matched samples extraction solution. Using other diluents may result in wrong results.

8. The solution and test card must be equilibrated to room temperature (20°C ~30°C) before used, otherwise the results may be incorrect.

9. Sensitivity maybe decrease if the sample did not test directly. Please test the sample as soon as possible.

10. Cross reactions maybe exist due to the N protein in SARS has a high homology with the new coronavirus (SARS-CoV-2). However, the interpretation of the results is not affected during seasons without SARS infection.

11. Analysis the possibility of false negative results:

1) Inappropriate sample collection, using other non-matching solution, sample transfer time is too long (more than half an hour), the volume of solution added when eluted the swab are too much, non- standardized elution operation, low virus titer in the sample, these may all lead to false negative results.

2) Mutations in viral genes may lead to changes in antigen epitope, leading to false negative results.

12. Analysis the possibility of false positive results:

1) Inappropriate sample collection, using other non-matching solutions, non-standardized elution operation, these may all lead to false positive results.

2) Cross-contamination of samples may lead to false positive results.

3) False negative result from nucleic acid.

13. Analysis the possibility of invalid result:

1) If the sample volume is not enough, the chromatography cannot be carried out successfully.

2) The test card would invalid if the package was broken. The packaging status must be carefully checked before use.

14. In different stages of infection, the coincidence rate with nucleic acid detection results may be different due to different viral load.

15. When sampling a nasopharyngeal swab, both nostrils need to be sampled with the same swab. If you only take it once, it may cause wrong results.

16. The results may be different due to the different viral load of the same person in different time periods.

Performance Characteristics

SARS-CoV-2 :

1.Positive specificity

The results of the tests were positive for the SARS-CoV-2 antigen positive reference samples.

2.Negative specificity

The results of the tests were negative for the SARS-CoV-2 antigen negative reference samples.

3.Limit of detection

The results of the tests were in line with the requirements of the minimum detection limit of SARS-CoV-2 antigen reference product.

4.Repeatability

The results were all positive and uniform in coloration after parallel testing for 10 times with reference samples of SARS-CoV-2 antigen.

Influenza A:

1.Positive specificity

The results of the tests were positive for the influenza A virus antigen positive reference samples.

2.Negative specificity

The results of the tests were negative for the influenza A virus antigen negative reference samples.

3.Limit of detection

The results of the tests were in line with the requirements of the minimum detection limit of influenza A virus antigen reference product.

4.Repeatability

The results were all positive and uniform in coloration after parallel testing for 10 times with reference samples of influenza A virus antigen.

Influenza B:

1.Positive specificity

The results of the tests were positive for the influenza B virus antigen positive reference samples.

2.Negative specificity

The results of the tests were negative for the influenza B virus antigen negative reference samples.

3.Limit of detection

The results of the tests were in line with the requirements of the minimum detection limit of influenza B virus antigen reference product.

4.Repeatability

The results were all positive and uniform in coloration after parallel testing for 10 times with reference samples of influenza B virus antigen.

Cross-reactivity

SARS-CoV-2 :

The results showed no cross reactivity with influenza a virus, influenza B virus, respiratory adenovirus, respiratory syncytial virus and mycoplasma pneumoniae.

Influenza A:

The results showed no cross reactivity with respiratory adenovirus, respiratory syncytial virus and mycoplasma pneumoniae.

Influenza B:

The results showed no cross reactivity with respiratory adenovirus, respiratory syncytial virus and mycoplasma pneumoniae.

Interfering

The test result of SARS-CoV-2 and Influenza A+B Antigen Rapid Test do not be interfered with the following drugs : zanamivir, ribavirin, oseltamivir, levofloxacin cefradine meropenem, tobramycin, oxymetazoline hydrochloride nasal spray, budesonide.

Precaution

1. The reagent is a disposable diagnostic reagent in vitro, which is only used for the detection of human nasopharyngeal swab, or oropharyngeal swab. The operation should be carried out strictly according to the instructions. Do not use expired and damaged products.

2. The strength of the quality control line does not mean the quality of the reagent, as long as its color is clear and visible, that means the reagent is effective.

3. The kit should be sealed and kept away from moisture. Reagents or samples stored at low temperature should be balanced to room temperature before they can be used.

4. Reagents should be used as soon as possible after removal from aluminum

foil bags, so as to avoid exposure to air for too long and affecting test results due to dampness.

5. Do not use samples that have been placed for too long or contaminated.

6. Please operate in accordance with the laboratory testing procedures for infectious diseases. Waste after use should be treated in accordance with infectious substances and should not be discarded at will.

Note: use clean pipettes or nozzles for each sample to avoid cross contamination.

7.Incorrect operation may affect the accuracy of the results, such as insufficient or excessive sample extract, insufficient sample mixing, insufficient sample volume, inaccurate detection time, etc.

8. Components in different batch should not be mixed.

9.If the sample swab is not rotated and squeezed in the sample extraction tube for 10 times, false negative results may occur. If the swab is put into the packaging bag after sample collection, false negative results may occur.

10. There should be appropriate biosafety assurance procedures for those substances containing and suspected sources of infection. The following are relevant considerations:

1) Handle samples and reagents with gloves;

2) Do not suck samples with your mouth;

3) Do not smoke, eat, drink, cosmetic or handle contact lenses while handling these items;

4) Disinfect the spilled sample or reagent with disinfectant;

5) Disinfect and treat all samples, reagents and potential pollutants in accordance with relevant local regulations;

Each component of the reagent remains stable until the expiry date under proper handling and storage conditions. Do not use the expired reagent kit.

MANUFACTURER / POST-SALE SERVICE UNIT

Weifang Kanghua Biotech Co., Ltd.

Add.: No.669 Yuehe Road, Economic Industrial Development Zone, Weifang, Shandong, China

EUROPEAN REPRESENTATIVE

CMC Medical Devices & Drugs S.L.

Add: C/ Horacio Lengo N° 18, CP 29006, Málaga, Spain



INSTRUCTIONS OF SYMBOL

	Consult instructions for use		Keep dry
	Temperature limit	LOT	Batch code
	For single use	IVD	In vitro diagnostic medical device
	Manufacturer		Date of manufacture
	Use-by date		Contains sufficient for <n> tests
	Keep away from sunlight		

IFU SARS-CoV-2/Influ A+B/ Antigen A/O